

# Catalyst Dx\* Chemistry Analyser

## QC/Maintenance Log

The Catalyst Dx\* analyser performs a self-check every time it is powered on, when a test is performed, and when it is initialised.

Year: \_\_\_\_\_

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Clean the internal components of the analyser.												
	Run VetTrol* quality control.												
	Clean the centrifuge.												
	Clean the waste drawer.												
Quarterly	Clean the fan filter.												
When Prompted	Clean the internal components of the analyser.												
As Needed	Clean the case.												
	Clean the screen.												
	Clean the waste drawer.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the Catalyst Dx Analyser, see the *Catalyst Dx Chemistry Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# Catalyst One\* Chemistry Analyser

## QC/Maintenance Log

The Catalyst One\* analyser performs a self-check every time it is powered on, when a test is performed, and when it is initialised.

Year: \_\_\_\_\_

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Clean the internal components of the analyser.												
	Run VetTrol* quality control.												
	Clean the waste drawer.												
When Prompted	Clean the internal components of the analyser.												
As Needed	Clean the case.												
	Clean the waste drawer.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the Catalyst One Analyser, see the *Catalyst One Chemistry Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# Coag Dx\* Analyser

## QC/Maintenance Log

The Coag Dx\* Analyser performs a self-check every time it is activated, when a test is performed, and when an electronic quality control is run daily.

Month: \_\_\_\_\_

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
As Needed	Clean the key pad, exterior case, and screen.																																

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the Coag Dx Analyser, see the *Coag Dx Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# LaserCyte\* Hematology Analyser

## QC/Maintenance Log

The LaserCyte analyser is a liquid-based system that needs power to maintain its cleansing cycle. If the LaserCyte analyser is powered off for more than 24 consecutive hours, please contact IDEXX Customer Support.

Quarter: \_\_\_\_\_

Frequency	Maintenance Required	Week 1	Week 2	Week 3	Week 4	Week 5	Week 6	Week 7	Week 8	Week 9	Week 10	Week 11	Week 12
Weekly	Clean the air filter.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the LaserCyte Analyser, see the *LaserCyte Haematology Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# ProCyte Dx\* Haematology Analyser

## QC/Maintenance Log

Month: \_\_\_\_\_

The analyser will perform a background check when it restarts after exiting Standby mode. It takes approximately 8 minutes to exit Standby mode.

By default, the ProCyte Dx analyser will set itself to Standby mode at 7:00 p.m. each day. To change the default setting, see the Operator's Guide or contact IDEXX Customer and Technical Support.

To power off the analyser, tap **Turn Off ProCyte Dx** on the ProCyte Dx Instruments screen and follow the on-screen instructions. When powering on the IPU and ProCyte Dx analyser, ensure that the IDEXX VetLab\* Station is powered back on first.

Clean the fan filter while the analyser is powered off for the weekly reboot.

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
As Needed	Change the reagent kit.																															
	Change the stain pack.																															
Weekly	Power off the analyser and IPU.																															

Year: \_\_\_\_\_

The Monthly Rinse procedure requires a 5% bleach solution of filtered, unscented bleach with distilled/deionised water. The ratio of bleach to water depends on the bleach concentration. (For example, Clorox\* concentrated bleach has a 7.5% concentration, so the solution should be 2 parts bleach with 1 part deionised/distilled water.)

Use appropriate lab handling techniques when using e-CHECK\* (XS) quality control.

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Clean the fan filter.												
	Perform the Monthly Rinse procedure (when prompted by the IDEXX VetLab Station).												
	Run e-CHECK (XS) L2 quality control.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the ProCyte Dx analyser, see the *ProCyte Dx Haematology Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# ProCyte One\* Haematology Analyser

## QC/Maintenance Log

Month: \_\_\_\_\_

The ProCyte One\* Haematology Analyser performs most maintenance automatically outside of business hours. It runs a self-cleaning sequence after each sample run and cycles fluid daily outside of normal business hours. A more robust weekly maintenance sequence runs on Sundays at 3:00 a.m., including system checks, rinses, and an automated IDEXX SmartQC\* control run. IDEXX will occasionally send software updates to the ProCyte One. No action is required on your part to install the updates.

**IMPORTANT:** The ProCyte One must remain powered on to receive software updates and to allow automated maintenance cycles to run. The analyser will go into low-power mode after 15 minutes of idle time and will reboot periodically as needed during weekly maintenance.

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
As Needed	Change the reagent.																																
As Needed	Change the sheath.																																
As Needed	Change IDEXX SmartQC*.																																

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the ProCyte One analyser, see the *ProCyte One Haematology Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# IDEXX SediVue Dx\* Urine Sediment Analyser

## QC/Maintenance Log

Month: \_\_\_\_\_

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
As Needed	Empty the waste bin.																															
Weekly	Reboot the analyser.																															

Year: \_\_\_\_\_

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Clean the analyser.												
	Clean the fan filter.												
	Run L1 and L2 quality control.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the SediVue Dx analyser, see the *IDEXX SediVue Dx Urine Sediment Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# SNAP Pro\* Analyser

## QC/Maintenance Log

Month: \_\_\_\_\_

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Weekly	Power the device off.																																
Weekly/As Needed	Clean the base, reflectors, and activator.																																

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the SNAP Pro analyser, see the *SNAP Pro Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).





# SNAPshot Dx\* Analyser

## QC/Maintenance Log

The SNAPshot Dx Analyser calibrates automatically. It performs a self-check at the beginning of each test as well as an automatic systems check each time the analyser is powered on.

Year: \_\_\_\_\_

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
As Needed	Clean the case.												
	Clean the SNAP* ports.												
	Clean the screen.												
	Clean the incubator wells.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the SNAPshot Dx Analyser, see the *SNAPshot Dx Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# IDEXX VetAutoread\* Haematology Analyser

## QC/Maintenance Log

The IDEXX VetAutoread\* analyser performs a self-check every time it is powered on and whenever a Cal Rod is run.

Month: \_\_\_\_\_

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Daily	Calibrate the analyser with a Cal Rod.																															

For more information about maintaining the VetAutoread Analyser, see the *VetAutoread Haematology Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# IDEXX VetLab\* Station

## QC/Maintenance Log

Quarter: \_\_\_\_\_

Frequency	Maintenance Required	Week 1	Week 2	Week 3	Week 4	Week 5	Week 6	Week 7	Week 8	Week 9	Week 10	Week 11	Week 12
Weekly	Back up data.												
	Reboot the computer.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the IDEXX VetLab Station, see the *IDEXX VetLab Station Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# IDEXX VetLab\* UA\* Analyser

## QC/Maintenance Log

Month: \_\_\_\_\_

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Daily	At end of day, clean the test strip tray† with water and, if necessary, 70% isopropyl alcohol.																															
Weekly	Calibrate the analyser every seven days.																															

†**IMPORTANT:** Be careful not to damage the grey reference pad during cleaning, and ensure that it is completely clean and dry before reading a test strip.

For more information about maintaining the IDEXX VetLab UA Analyser, see the *IDEXX VetLab UA Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# VetStat\* Electrolyte and Blood Gas Analyser

## QC/Maintenance Log

The VetStat analyser performs a self-check every time it is activated, when a test is performed, and when an electronic quality control is run daily.

Month: \_\_\_\_\_

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Daily	Run a standard reference cassette (SRC), alternate levels.																															
Weekly	Clean the sample measurement chamber.																															

Year: \_\_\_\_\_

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Run OPTI* CHECK quality control (alternate levels).												
Quarterly	Run a haemoglobin calibration cassette (HbCC) and reinstall the standard reference cassettes (i.e., rescan SRC setup barcodes).												
Annually	Replace peristaltic pump cartridge.												
	Change gas I/O port.												

For more information about maintaining the VetStat Analyser, see the *VetStat Electrolyte and Blood Gas Analyser Operator's Guide* at [idexx.com/product-info](http://idexx.com/product-info).



# VetTest\* Chemistry Analyser

## QC/Maintenance Log

The VetTest analyser performs a self-check every time it is activated, when a test is performed, and when an electronic quality control is run daily.

Year: \_\_\_\_\_

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Run VetTrol* quality control with the quality control panel.												
	Clean pipettor actuator collar with an alcohol swab.												

### Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the VetTest Analyser, see the *VetTest Chemistry Analyser Operator's Guide*.